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Materasso Antidecubito Oba 1T.

Codice	Misure	Portata Min	Portata Max
05.05.OBA1TTS	200x90x15 cm	40 Kg	150 Kg
05.05.OBA1TTS120	200x120x15 cm	40 Kg	150 Kg

**Destinazione d'uso**

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**Use**

Device intended for the prevention of bedsores in subjects who must maintain a lying position for prolonged periods of time, geriatric patients, subjects with poor mobility and subjects with sensitive and red skin.

Recommended for pediatric use and for patients weighing less than 50 Kg. Intended for use in hospital settings, pediatric care facilities, retirement homes and nursing homes. The device is classified with anti-decubitus grade 3, according to the EPUAP scale.

**Description**

The OBA 1T mattress is made of EvoPore foam with differentiated density and processing.

The mattress guarantees effective prevention of the formation of pressure ulcers and ensures pleasant comfort by providing ideal support thanks to the foam that adapts to the morphology of the patient's body.

The reinforced edges guarantee good stability, reducing the risk of slipping, making it easier for the person to lift, helping them get in and out of bed and excellently support the patient's lateral transfers. The mattress is equipped with two symmetrical surfaces with the same workmanship, which allow the mattress to be turned over cyclically, significantly extending its duration of use.

The EvoPore foam maintains its very high stability and high support power even during the night and stands out for its excellent resistance and durability over time. The lining is made of elastic, soft and resistant TRICOT SANILON cotton with a zip closure on three sides.

Materials	
Support Layer	EvoPore Foam with a Density of 35 Kg/m <sup>3</sup>
Intermediate Layer	EvoPore Foam with a Density of 30 Kg/m <sup>3</sup>
Base	EvoPore Foam with a Density of 35 Kg/m <sup>3</sup>
Borders	EvoPore Foam with a Density of 40 Kg/m <sup>3</sup>
Cover	Tricot Sanilon
Cover Maintenance and Cleaning	
	Water Temperature not Above 40°C , Delicate
	Do not Bleach
	Maximum Temperature 150°C or 300°F
	Dryclean with any Solvent

Tumble Dry Medium Heat

### Mattress Maintenance and Cleaning

Possible cleaning with a damp cloth and neutral detergent

Before use dry the mattress carefully without exposing it to direct sunlight

### Storage and Disposal

Single-user Product

Keep away from heat, sparks, open flame or other sources of ignition

Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

### Package and Nr.Items

Weight of the mattress

10 Kg

Box

1 Piece

### Warranty

European

24 months for manufacturing defects with the exception of consumable parts

Oba Warranty

10 years as reported in the user manual

### Product References

Manufacturer

OBA AG

Authorized Distributor

Demarta Virginio S.r.l.

Directive MDR 2017/745

Class I

ISO Code DM 332/99

03.33.06.006

ISO Code DPCM 12.01.2017

04.33.06.006

CND Code

V080702

RDM Code

2545057/R