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Fascia Addominale Beauty Line.

Destinazione d'uso

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Waist Cuff Beauty Line.

Use

The waist cuffs of the Beauty Line series are combinable with 07.01.HF3006 and 07.01.HF3008 pumps. They are devices for lymphatic drainage; a treatment against water retention. The progressive pressure stimulate the flow of body fluids, facilitating their elimination. The pressure therapy is a medical device, do not use without a physician opinion in case of health problems.

Features

Code	07.01.WAIST12037
------	------------------

Material	Nylon • PU
----------	------------

Warnings

Pressure therapy is contraindicated for people suffering from the following pathologies: cardiopathy, suspected or ascertained thrombosis, arteriosclerosis or ischemia, dermatitis, gangrene, local infections

Check the correct assembly of the device with particular attention to the connection of the tubes between the pump and the waist cuff

Maintenance and Cleaning

Clean the waist cuff with a neutral soap; the drying must not take place in direct sunlight

Storage and Disposal

Single-user Product

Keep away from heat, sparks, open flame or other sources of ignition

Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

Package and Nr. Items

Box	1 Piece
-----	---------

Warranty

European	24 months for manufacturing defects with the exception of consumable parts
----------	--

Product References

Manufacturer	Demarta Virginio S.r.l.
--------------	-------------------------

Directive MDR 2017/745	Class I
------------------------	---------

CND Code	Z120607
----------	---------

RDM Code	1745140/R
----------	-----------