



### Nome Prodotto

Posizionatore per Braccio.

| Codice      | Misure (cm)   |
|-------------|---------------|
| 05.00.AP071 | 40 x 12.5 x 2 |
| 05.00.AP072 | 60 x 12.5 x 2 |

### Categoria

Posizionatore in Gel per Sala Operatoria.

### Descrizione

Il POSIZIONATORE PER BRACCIO protegge e stabilizza gli arti superiori durante gli interventi chirurgici.

### Materiali

Il POSIZIONATORE PER BRACCIO è realizzato in gel amorfo di poliuretano con densità di 0,99 gr/cm<sup>3</sup> rivestito da un film di poliuretano.

### Manutenzione e Pulizia

- I posizionatori in gel posso essere puliti facilmente con qualsiasi detergente e/o disinfettante per la sala operatoria.
- Evitare di immergere per lungo tempo i posizionatori in gel nel detergente e/o disinfettante.

### Conservazione e Smaltimento

- Mantenere, se possibile, i posizionatori in una superficie in piano.
- Evitare di appoggiare sui posizionatori oggetti pesanti o taglienti.
- Mantenere la temperatura del gel tra i 50° e i -18°.
- Conservare al riparo da fonti di luce, umidità e calore.
- Smaltire il dispositivo rispettando le normative vigenti in materia di tutela ambientale e raccolta differenziata.

### Confezione

1 Pezzo.

### Garanzia

24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili.

### Codice ISO

04.33.09

### Codice CND


Y033399

### Codice RDM

90736/R

### Classe Dispositivo Medico

I conforme alla Direttiva 93/42/CEE



Product Name  
Armboard Pad.



| Code        | Dimensions (cm) |
|-------------|-----------------|
| 05.00.AP071 | 40 x 12.5 x 2   |
| 05.00.AP072 | 60 x 12.5 x 2   |

Category  
Gel Positioner for Operating Table.

Description  
It can be used as arm support during all types of surgery.

Materials  
It is made of an amorphous polyurethane gel with a density of 0.99 gr/cm<sup>3</sup>.

Maintenance and Cleaning

- The gel positioners can be easily cleaned with any cleaning agent and / or disinfectant for the operating room.
- Avoid immersing the gel positioners in the detergent and / or disinfectant for a long time.

Storage and Disposal

- Keep the positioners, if possible, on a level surface.
- Avoid placing heavy or sharp objects on the positioners.
- Keep the gel temperature between 50 ° and -18 °.
- Keep away from sources of light, moisture and heat.
- Dispose of the device in compliance with current regulations on environmental protection and separate waste collection.

Package  
1 Piece.

Warranty  
24 months for manufacturing defects with the exception of consumable parts.

ISO Code  
04.33.09

CND Code  
Y033399

RDM Code  
90736/R

Medical Device Class  
I complies with Directive 93/42/CEE

