

**Nome Prodotto**

Carrozzina da Transit Viky Transit Plus.

**Destinazione d'uso**

Le carrozzine sono un mezzo di trasporto indispensabile per le persone disabili con assente o ridotta capacità di deambulazione.

Lo scopo è quello di facilitare gli spostamenti e di salvaguardare il benessere fisico della persona, evitando la comparsa di complicanze correlate alla posizione seduta obbligatoria (piaghe da decubito, dolori posturali cronici...).



Codice	Taglia	Lunghezza	Larghezza	Altezza	Larghezza Piegata	Larghezza Seduta	Profondità Seduta	Altezza Schienale	Altezza Seduta da Terra	Portata Massima (Kg)	Dimensioni Ruote Anteriori	Dimensioni Ruote Posteriori	Peso Netto con Poggiapiedi	Dimensioni Scatola
Misure in Cm														
01.04.320D41FD	41	69	53.5	92	22.5	41	40	44	50	120	20 (8")	30 (12")	15	72x25x94
01.04.320D43FD	43	69	55.5	92	23	43	40	44	50	120	20 (8")	30 (12")	15.1	72x25x94
01.04.320D46FD	46	69	57.5	92	23.5	46	40	44	50	120	20 (8")	30 (12")	15.2	72x25x94
01.04.320D48FD	48	69	59.5	92	24	48	40	44	50	120	20 (8")	30 (12")	15.3	72x25x94
01.04.320D50FD	50	69	61.5	92	25	50	40	44	50	120	20 (8")	30 (12")	15.4	72x25x94

**Caratteristiche**

Telaio	Acciaio Verniciato a Polvere
Crociera	Doppia
Colore	Grigio Canna di Fucile (Opaco)
Ruote Anteriori	20 x 2,5 cm (8" x 1") PVC Grigio
Ruote Posteriori	30 cm (12") PU Grigio
Braccioli	Ribaltabili, Modello a Scrivania con Poggiabraccio in PVC
Pedane	Estraibili e Regolabili con Supporto per Gambe
Poggiapiedi	PVC Nero
Sedile e Schienale	Nylon Morbido Imbottito, lavabile, disinfettabile e ignifugo

**Manutenzione e Pulizia**

Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica

**Conservazione e Smaltimento**

Si consiglia l'uso ad un solo utente

Conservare al riparo da fiamme libere e fonti di calore

Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata

**Confezione e Nr. Articoli**

Scatola 1 Pezzo

**Garanzia**

Europea 24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili

**Riferimenti del Prodotto**

Fabbricante Demarta Virginio S.r.l.

Classe di Dispositivo MDR 2017/745 Classe I

Codice ISO DM332/99 12.21.06.045 / 12.24.06.175

Codice ISO DPCM 12.01.2017 12.22.18.003

Codice CND Y122106

Codice RDM

2059397/R















### Product Name

Transport Wheelchairs Viky Transit Plus.

### Use

The wheelchairs are an essential way of transport for disabled people with none or limited walking ability. The purpose is to facilitate movements and to safeguard the physical well-being of the person, avoiding the appearance of complications related to the sitting position (bed sores, chronic postural pains ...).



Code	Size	Length	Width	Height	Width Folded	Seat Width	Seat Depth	Back Height	Floor to Seat Height	Weight Capacity	Front Wheel	Rear Wheel	Weight with Footrest	Carton Size
Measures in Cm														
01.04.320D41FD	41	69	53.5	92	22.5	41	40	44	50	120	20 (8")	30 (12")	15	72x25x94
01.04.320D43FD	43	69	55.5	92	23	43	40	44	50	120	20 (8")	30 (12")	15.1	72x25x94
01.04.320D46FD	46	69	57.5	92	23.5	46	40	44	50	120	20 (8")	30 (12")	15.2	72x25x94
01.04.320D48FD	48	69	59.5	92	24	48	40	44	50	120	20 (8")	30 (12")	15.3	72x25x94

### Features

Frame	Powder Coating Steel Frame
Cross Bar	Double Cross Brace
Color	Smokey Grey (Dull Polish)
Front Wheels	20 x 2,5 cm (8" x 1") PVC Grey Tyre
Rear Wheels	30 cm (12"), Spoke Wheels, PU Grey Tyre
Armrest	Flip Back Desk Model Armrest with PVC Armpad
Footrest	Swing Away and Length Adjustable with Leg Support
Footplate	Plastic Black
Seat and Backrest	Soft and Padded Nylon, washable, disinfectable and fire-retardant

### Maintenance and Cleaning

Disinfection is possible with 5% diluted disinfectants even with alcohol base

### Storage and Disposal

Single-user Product

Keep away from heat, sparks, open flame or other sources of ignition

Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

### Package and Nr. Items

Box 1 Piece

### Warranty

European 24 months for manufacturing defects with the exception of consumable parts

## Product References

Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Class I
ISO Code DM332/99	12.21.06.045 / 12.24.06.175
ISO Code DPCM 12.01.2017	12.22.18.003
CND Code	Y122106
RDM Code	2059397/R