

**Nome Prodotto**

Fodera in Nylon per Minislide in cotone.

Codice	Misure (cm)	Abbinamento
03.04.2051	56 x 63	03.04.2030

**Destinazione d'uso**

I Sistemi di Trasferimento sono destinati alla facilitazione delle operazioni di movimentazione manuale dei pazienti attraverso manovre di spinta, trazione, scivolamento e scorrimento. Sono stati progettati per l'uso in ospedali, case di cura e strutture sanitarie/assistenziali e per l'assistenza domiciliare.

Descrizione

La FODERA IN NYLON è pensata per poter rendere impermeabile, facile da pulire e per evitare un'eccessiva usura dei seguenti sacchi di trasferimento:

- Minislide in cotone, codice: 03.04.2030

Materiali	
Fodera	Nylon
Manutenzione e Pulizia	
	Lavaggio Max. 60°C
Disinfezione possibile con LT-CIN (45% alcool isopropilico)	
Conservazione e Smaltimento	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo MDR 2017/745	Classe I
Codice CND	Z1209080
Codice RDM	2504980

Product Name

Nylon Cover for cotton Minislide.

Code	Dimensions (cm)	Pairings
03.04.2051	56 x 63	03.04.2030



Use

The Transfer Systems are intended for facilitating the manual handling operations of patients through thrust, traction, sliding maneuvers. They have been designed for use in hospitals, nursing homes and healthcare facilities and for home care.

Description

The NYLON COVER is designed to make the following Sliding Mats waterproof and easy to clean:

- Cotton Minislide, Code: 03.04.2030

Materials	
Cover	Nylon
Maintenance and Cleaning	
	Water Temperature not Above 60°C
Disinfection Possible with LT-CIN (45% isopropyl alcohol)	
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Class I
CND Code	Z1209080
RDM Code	2504980