

Nome Prodotto

Teli di Scorrimento Sheet.



Codice	Misure (cm)	Nr maniglie	Codice	Misure (cm)	Nr maniglie	Codice	Misure (cm)	Nr maniglie
03.04.11429	195 x 70	31	03.04.11438	200 x 105	33	03.04.11450	95 x 75	17
03.04.11430	190 x 75	31	03.04.11439	130 x 80	23	03.04.11451	110 x 90	20
03.04.11431	170 x 72	27	03.04.11440	130 x 100	24	03.04.11452	70 x 50	13
03.04.11432	165 x 70	27	03.04.11441	170 x 90	28	03.04.11453	65 x 55	12
03.04.11433	150 x 70	25	03.04.11442	180 x 70	29	03.04.11454	110 x 75	19
03.04.11434	130 x 70	23	03.04.11443	175 x 110	31	03.04.11455	200 x 84	33
03.04.11435	110 x 80	19	03.04.11444	200 x 140	36	03.04.11456	170 x 60	26
03.04.11436	90 x 50	15	03.04.11448	185 x 90	31	03.04.11457	60 x 70	14
03.04.11437	200 x 70	31	03.04.11449	70 x 70	15	03.04.11458	40 x 50	10
03.04.11459	190 x 65	31	03.04.11460	80 x 100	17	03.04.11461	140 x 68	24





Destinazione d'uso

I Sistemi di Trasferimento sono destinati alla facilitazione delle operazioni di movimentazione manuale dei pazienti attraverso manovre di spinta, trazione, scivolamento e scorrimento. Sono stati progettati per l'uso in ospedali, case di cura e strutture sanitarie/assistenziali e per l'assistenza domiciliare.

Descrizione

Il telo di scorrimento SHEET è ideato per facilitare il posizionamento dei pazienti allettati. Si consiglia l'utilizzo contemporaneo di due SHEET: uno più grande a contatto del materasso e uno più piccolo a contatto del paziente, in modo tale che possano scorrere l'uno sull'altro. Tramite le numerose maniglie poste su tre lati, gli operatori possono muovere il paziente in tutte le direzioni con estrema facilità.

Portata Massima: 300 Kg.

Materiali	
Telo	Nylon Ultra Low Friction
Manutenzione e Pulizia	
	Lavaggio Max. 60°C
	Non Stirare
	Non Asciugatura a Tamburo
	Non è Possibile il Lavaggio a Secco
Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica	
Conservazione e Smaltimento	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo UE 2017/745	Classe I
Codice ISO DPCM 12.01.2017	12.31.03.009 / 12.31.03.012
Codice CND	Z12019010
Codice RDM	2468245/R

Product Name

Sliding Sheets • Sheet.







Code	Dim. (cm)	Nr handles	Code	Dim. (cm)	Nr handles	Code	Dim. (cm)	Nr handles
03.04.11429	195 x 70	31	03.04.11438	200 x 105	33	03.04.11450	95 x 75	17
03.04.11430	190 x 75	31	03.04.11439	130 x 80	23	03.04.11451	110 x 90	20
03.04.11431	170 x 72	27	03.04.11440	130 x 100	24	03.04.11452	70 x 50	13
03.04.11432	165 x 70	27	03.04.11441	170 x 90	28	03.04.11453	65 x 55	12
03.04.11433	150 x 70	25	03.04.11442	180 x 70	29	03.04.11454	110 x 75	19
03.04.11434	130 x 70	23	03.04.11443	175 x 110	31	03.04.11455	200 x 84	33
03.04.11435	110 x 80	19	03.04.11444	200 x 140	36	03.04.11456	170 x 60	26
03.04.11436	90 x 50	15	03.04.11448	185 x 90	31	03.04.11457	60 x 70	14
03.04.11437	200 x 70	31	03.04.11449	70 x 70	15	03.04.11458	40 x 50	10
03.04.11459	190 x 65	31	03.04.11460	80 x 100	17	03.04.11461	140 x 68	24

Use

The Transfer Systems are intended for facilitating the manual handling operations of patients through thrust, traction, sliding maneuvers. They have been designed for use in hospitals, nursing homes and healthcare facilities and for home care.

Description

The SLIDING SHEET is designed to facilitate the positioning of patients. It is recommended to use two sheets at the same time: a larger one in contact with the mattress and a smaller one in contact with the patient, so that they can slide over one another. Thanks to the numerous handles placed on three sides, the operators can move the patient in all directions with extreme ease. Maximum Capacity: 300 Kg.

Materials	
Sheet	Nylon Ultra Low Friction
Maintenance and Cleaning	
	Water Temperature not Above 60°C or 140°F
	Do not Iron
	Do not Tumble Dry
	Do not Dryclean
Disinfection is possible with 5% diluted disinfectants even with alcohol base	
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.r.l.
Directive UE 2017/745	Class I
ISO Code DPCM 12.01.2017	12.31.03.009 / 12.31.03.012
CND Code	Z12019010
RDM Code	2468245/R