

Nome Prodotto















Carrozzina ad Autospinta Viky Double Cross.

Destinazione d'uso

Le carrozzine sono un mezzo di trasporto indispensabile per le persone disabili con assente o ridotta capacità di deambulazione.

Lo scopo è quello di facilitare gli spostamenti autonomi e di salvaguardare il benessere fisico della persona, evitando la comparsa di complicanze correlate alla posizione seduta obbligatoria (piaghe da decubito, dolori posturali cronici...).



Codice	Taglia	Lunghezza	Larghezza	Altezza	Larghezza Piegata	Larghezza Seduta	Profondità Seduta	Altezza Schienale	Altezza Seduta da Terra	Portata Massima (Kg)	Dimensioni Ruote Anteriori	Dimensioni Ruote Posteriori	Peso Netto con Poggiapiedi	Dimensioni Scatola
Misure in Cm														
01.04.310D41FD	41	80	56	90.5	26	41	40	44	50	120	20 (8")	60 (24")	18.7	82x27x93
01.04.310D43FD	43	80	58.5	90.5	26.5	43	40	44	50	120	20 (8")	60 (24")	18.8	82x27x93
01.04.310D46FD	46	80	61	90.5	27	46	40	44	50	120	20 (8")	60 (24")	18.9	82x27x93
01.04.310D48FD	48	80	64	90.5	27.5	48	40	44	50	120	20 (8")	60 (24")	19	82x27x93
01.04.310D50FD	50	80	65	90.5	28	50	40	44	50	120	20 (8")	60 (24")	20	82x27x93

Caratteristiche	
Telaio	Acciaio Verniciato a Polvere
Crociera	Doppia
Colore	Grigio Canna di Fucile (Opaco)
Ruote Anteriori	20 x 2,5 cm (8" x 1") PVC Grigie
Ruote Posteriori	60 cm (24") PU Grigie
Anello Corrimano	Acciaio Verniciato a Polvere Nero
Braccioli	Ribaltabili, Modello a Scrivania con Poggiabracchio in PVC
Pedane	Estraibili e Regolabili con Supporto per Gambe
Poggiapiedi	PVC Nero
Sedile e Schienale	Nylon Morbido Imbottito, lavabile, disinfettabile e ignifugo
Manutenzione e Pulizia	
Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica	
Conservazione e Smaltimento	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo MDR 2017/745	Classe I
Codice ISO DM332/99	12.21.06.033 / 12.21.06.039 / 12.24.06.175
Codice ISO DPCM 12.01.2017	12.22.03.003
Codice CND	Y122106
Codice RDM	2594008/R















Product Name

Manual Wheelchairs Viky Double Cross.

Use

The wheelchairs are an essential way of transport for disabled people with none or limited walking ability. The purpose is to facilitate autonomous movements and to safeguard the physical well-being of the person, avoiding the appearance of complications related to the sitting position (bed sores, chronic postural pains ...).



Code	Size	Length	Width	Height	Width Folded	Seat Width	Seat Depth	Back Height	Floor to Seat Height	Weight Capacity	Front Wheel	Rear Wheel	Weight with Footrest	Carton Size
Measures in Cm														
01.04.310D41FD	41	80	56	90.5	26	41	40	44	50	120	20 (8")	60 (24")	18.7	82x27x93
01.04.310D43FD	43	80	58.5	90.5	26.5	43	40	44	50	120	20 (8")	60 (24")	18.8	82x27x93
01.04.310D46FD	46	80	61	90.5	27	46	40	44	50	120	20 (8")	60 (24")	18.9	82x27x93
01.04.310D48FD	48	80	64	90.5	27.5	48	40	44	50	120	20 (8")	60 (24")	19	82x27x93
01.04.310D50FD	50	80	65	90.5	28	50	40	44	50	120	20 (8")	60 (24")	20	82x27x93

Features	
Frame	Powder Coating Steel Frame
Cross Bar	Double Cross Brace
Color	Smokey Grey (Dull Polish)
Front Wheels	20 x 2,5 cm (8" x 1") PVC Grey Tyre.
Rear Wheels	60 cm (24"), Spoke Wheels, PU Grey Tyre
Handrim	Powder Coating Black Steel
Armrest	Flip Back Desk Model Armrest with PVC Armpad
Footrest	Swing Away and Length Adjustable with Leg Support
Footplate	Plastic Black
Seat and Backrest	Soft and Padded Nylon, washable, disinfectable and fire-retardant
Maintenance and Cleaning	
Disinfection is possible with 5% diluted disinfectants even with alcohol base	
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr. Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Class I
ISO Code DM332/99	12.21.06.033 / 12.21.06.039 / 12.24.06.175
ISO Code DPCM 12.01.2017	12.22.03.003
CND Code	Y122106
RDM Code	2594008/R

