


**Nome Prodotto**

Cuscino antidecubito in Gel ed Espanso Viscoelastico Temperfoam.

Codice	Misure (cm)	Portata (kg)
05.01.H3720-38	38 x 38 x 7	120
05.01.H3720-40	40 x 40 x 7	120
05.01.H3720-42	42 x 42 x 7	120
05.01.H3720-45	45 x 45 x 7	120
05.01.H3720-50	50 x 45 x 7	120

**Destinazione d'uso**

Dispositivo destinato alla prevenzione delle piaghe da decubito in soggetti che devono mantenere la posizione seduta per prolungati periodi di tempo.

DescrizioneIl cuscino TEMPERFOAM è composto da una base in espanso con densità di 30 Kg/m³ e da una superficie in espanso viscoelastico con densità di 60 Kg/m³. La fodera esterna è impermeabile e traspirante con una base antiscivolo e cerniera su 3 lati; i laccetti alle estremità permettono l'aggancio alla carrozzina.


Materiali	
Parte Inferiore Cuscino	Espanso con densità di 30 Kg/m ³
Parte Superiore Cuscino	Espanso viscoelastico con densità di 60 Kg/m ³
Fodera	25% Poliestere / 75% Poliuretano
Manutenzione e Pulizia Cuscino	
Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica	
Manutenzione e Pulizia Fodera	
	Lavaggio Max 95° C
	Non Stirare
	Non Candeggiare
	Non Utilizzare l'Asciugatrice
Conservazione e Smaltimento	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo MDR 2017/745	Classe I
Codice ISO DM332/99	03.33.03.009
Codice CND	V08030102
Codice RDM	2598408/R



Product Name

Viscoelastic Gel and Foam Cushion Temperfoam.

Code	Dimensions (cm)	Capacity (kg)
05.01.H3720-38	38 x 38 x 7	120
05.01.H3720-40	40 x 40 x 7	120
05.01.H3720-42	42 x 42 x 7	120
05.01.H3720-45	45 x 45 x 7	120
05.01.H3720-50	50 x 45 x 7	120







Use

Device designed for the prevention of bedsores in subjects who have to maintain the sitting position for a long time.

Description

The TEMPERFOAM cushion is composed of 30 Kg/m3 density foam base and 60 Kg/m3 density memory foam contact surface. The external cover is waterproof breathable with an anti-slip base and zip on three sides; 4 laces on the side allows to tie the cushion to wheelchairs and chairs.

Materials	
Cushion Lower Part	30 Kg/m3 density foam
Cushion Upper Part	60 Kg/m3 density memory foam
Cover	25% Polyester / 75% Polyurethane
Cushion Maintenance and Cleaning	
Disinfection is possible with 5% diluted disinfectants even with alcohol base	
Cover Maintenance and Cleaning	
	Water temperature not above 95° C
	Do not iron
	Do not bleach
	Do not dry in dryer
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Classe I
ISO Code DM332/99	03.33.03.009
CND Code	V08030102
RDM Code	2598408/R