

Nome Prodotto

Imbracature Deluxe 200 senza Testiera in Maglia.

Codice	Misure	Portata (Kg)
01.01.2022.200	S	200
01.01.2023.200	M	200
01.01.2024.200	L	200
01.01.2025.200	XL	200




**Destinazione d'uso**

Le Imbracature vengono associate ad un sollevatore della Linea SCANMEDICALEUROPE per la movimentazione dei pazienti.

Descrizione

L'IMBRACATURA DELUXE 200 SENZA TESTIERA, grazie al suo design, consente il trasferimento del paziente da seduto a sdraiato e viceversa.

Per garantire la massima sicurezza, le imbracature dovrebbero essere sostituite ogni 3 anni.

Materiali	
Imbracatura	Maglia
Manutenzione e Pulizia	
	Lavaggio Max. 40°C
	Non Stirare
	Non Candeggiare
	Non Asciugare a Tamburo
Conservazione e Smaltimento	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo MDR 2017/745	Classe I
Codice ISO DPCM 12.01.2017	12.36.21.003
Codice CND	V0805030102
Codice RDM	253031 1/R

Product Name

Mesh Deluxe 200 without Head Support.

Code	Size	Capacity (Kg)
01.01.2022.200	S	200
01.01.2023.200	M	200
01.01.2024.200	L	200
01.01.2025.200	XL	200







Use

The slings are associated with a lift of the SCANMEDICALEUROPE Line for patient handling.

Description

Its design allows transfers of patients who are in a sitting or lying position.
To ensure maximum safety, the sling should be replaced every 3 years.

Materials	
Sling	Mesh
Maintenance and Cleaning	
	Water Temperature not Above 40°C or 105°F
	Do not Iron
	Do not Bleach
	Do not Tumble Dry
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Class I
ISO Code DPCM 12.01.2017	12.36.21.003
CND Code	V0805030102
RDM Code	253031 1/R