

Nome Prodotto

Carrozzina da Transito Viky Transit Light.

Destinazione d'uso

Le carrozzine sono un mezzo di trasporto indispensabile per le persone disabili con assente o ridotta capacità di deambulazione.

Lo scopo è quello di facilitare gli spostamenti e di salvaguardare il benessere fisico della persona, evitando la comparsa di complicanze correlate alla posizione seduta obbligatoria (piaghe da decubito, dolori posturali cronici...).



Codice	Taglia	Lunghezza	Larghezza	Altezza	Larghezza Piegata	Larghezza Seduta	Profondità Seduta	Altezza Schienale	Altezza Seduta da Terra	Portata Massima (Kg)	Dimensioni Ruote Anteriori	Dimensioni Ruote Posteriori	Peso Netto con Appoggiatesta	Dimensioni Scatola
Misure in Cm														
01.04.640DQR41FD	41	97	56	92	26	41	41	43	49	120	20 (8")	30 (12")	12	74x27x94
01.04.640DQR43FD	43	97	58	92	26	43	41	43	49	120	20 (8")	30 (12")	12	74x27x94
01.04.640DQR46FD	46	97	61	92	26	46	41	43	49	120	20 (8")	30 (12")	12	74x27x94
01.04.640DQR48FD	48	97	63	92	26	48	41	43	49	120	20 (8")	30 (12")	12	74x27x94

*Modello 46 cm e 48 cm Iso: 12.24.06.175

Caratteristiche

Telaio	Alluminio Verniciato a Polvere (Iso: 12.24.03.121)
Crociera	Doppia (Iso: 12.24.03.118)
Colore	Grigio Canna di Fucile (Opaco)
Ruote Anteriori	20 cm (8") Piroettanti
Ruote Posteriori	30 cm (12")
Braccioli	Modello a Scrivania (Iso: 12.24.06.172), Ribaltabili e Regolabili in Altezza e Profondità (Iso: 12.24.06.160)
Pedane	Estraibili e Regolabili con Supporto per Gambe
Poggiatesta	PVC Nero
Sedile e Schienale	Nylon Morbido Imbottito, lavabile, disinfettabile e ignifugo
Ruote Anti-ribaltamento	

Manutenzione e Pulizia

Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica

Conservazione e Smaltimento

Si consiglia l'uso ad un solo utente

Conservare al riparo da fiamme libere e fonti di calore

Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata

Confezione e Nr. Articoli

Scatola 1 Pezzo

Garanzia

Europea 24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili

Riferimenti del Prodotto

Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo MDR 2017/745	Classe I
Codice ISO DM332/99	12.21.06.045 / 12.24.06.175
Codice ISO DPCM 12.01.2017	12.22.03.006 / 12.22.18.003
Codice CND	Y122106

Codice RDM

2459313/R















Product Name

Transport Wheelchairs Viky Transit Light.

Use

The wheelchairs are an essential way of transport for disabled people with none or limited walking ability. The purpose is to facilitate movements and to safeguard the physical well-being of the person, avoiding the appearance of complications related to the sitting position (bed sores, chronic postural pains ...).



Code	Size	Length	Width	Height	Width Folded	Seat Width	Seat Depth	Back Height	Floor to Seat Height	Weight Capacity	Front Wheel	Rear Wheel	Weight with Footrest	Carton Size
Measures in Cm														
01.04.640DQR41FD	41	97	56	92	26	41	41	43	49	120	20 (8")	30 (12")	12	74x27x94
01.04.640DQR43FD	43	97	58	92	26	43	41	43	49	120	20 (8")	30 (12")	12	74x27x94
01.04.640DQR46FD	46	97	61	92	26	46	41	43	49	120	20 (8")	30 (12")	12	74x27x94
01.04.640DQR48FD	48	97	63	92	26	48	41	43	49	120	20 (8")	30 (12")	12	74x27x94

Features

Frame	Powder Coating Aluminum Frame
Cross Bar	Double Cross Brace
Color	Smokey Grey (Dull Polish)
Front Wheels	20 cm (8"), Swiveling
Rear Wheels	30 cm (12")
Armrest	Height-adjustable Flip Back Slope Model Armrest
Footrest	Swing Away and Length Adjustable with Leg Support
Footplate	Plastic Black
Seat and Backrest	Soft and Padded Nylon, washable, disinfectable and fire-retardant
Anti Tipper Wheels	

Maintenance and Cleaning

Disinfection is possible with 5% diluted disinfectants even with alcohol base

Storage and Disposal

Single-user Product

Keep away from heat, sparks, open flame or other sources of ignition

Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

Package and Nr. Items

Box 1 Piece

Warranty

European 24 months for manufacturing defects with the exception of consumable parts

Product References

Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Class I
ISO Code DM332/99	12.21.06.045 / 12.24.06.175
ISO Code DPCM 12.01.2017	12.22.03.006 / 12.22.18.003
CND Code	Y122106
RDM Code	2459313/R