

Nome Prodotto  
Rollator Compacto in alluminio.

Codice	Colore
01.06.401A	Argento
01.06.401GO	Grigio Opaco
01.06.401B	Blu Notte



#### Destinazione d'uso

I deambulatori sono indicati nella fase di rieducazione alla deambulazione in soggetti anziani o con limitazioni funzionali di varia natura. Offrono stabilità, facilità di utilizzo e migliorano la sensazione di equilibrio dell'utilizzatore. Quest'ausilio è concepito per l'uso di un solo utente, in ambienti interni o esterni in base al modello.

	(Altezza sedile livello 1) *	(Altezza sedile livello 2) *	(Altezza sedile livello 3) *	(Altezza sedile livello 4) *
Altezza delle manopole (cm)	77 / 90 (in 8 posizioni)	79 / 92 (in 8 posizioni)	81 / 94 (in 8 posizioni)	83,5 / 96,5 (in 8 posizioni)
Altezza del sedile (cm)	50	52	54	56,5
Dimensioni sedile (cm)	36 x 31	36 x 31	36 x 31	36 x 31
Lunghezza telaio spiegato (cm)	68,5	71	74	76,5
Larghezza totale (cm)	62	62	62	62
Distanza le manopole (cm)	46	46	46	46
Dimensioni della ruota (cm)	20 (8")	20 (8")	20 (8")	20 (8")
Dimensioni borsa (cm)	33,5 x 21 x 15	33,5 x 21 x 15	33,5 x 21 x 15	33,5 x 21 x 15
Portata massima (kg)	110 kg	110 kg	110 kg	110 kg
Peso netto (kg)	7,2	7,2	7,2	7,2
Peso lordo (kg)	7,8	7,8	7,8	7,8
Dimensioni cartone (cm)	54,5 x 16 x 49,5	54,5 x 16 x 49,5	54,5 x 16 x 49,5	54,5 x 16 x 49,5

#### Descrizione

- Struttura in alluminio pieghevole.
- Seduta regolabile in 4 altezze \*.
- Seduta e schienale imbottiti con supporto laterale porta-bastone.
- Quattro ruote del diametro di 20 cm (ruote frontali piroettanti, due posteriori fisse con freno).
- Cestino porta oggetti incluso.
- Manopole ergonomiche e regolabili in altezza.
- Schienale imbottito regolabile in inclinazione.

Materiali	
Telaio	Alluminio
Impugnatura	Polipropilene
Ruote	PU
Seduta	PVC Imbottito
Manutenzione e Pulizia Telaio	
Usare un panno umido con del sapone neutro	
Non utilizzare macchine di lavaggio a getto di acqua e/o vapore, oggetti abrasivi e spazzole metalliche	
Evitare assolutamente l'utilizzo di prodotti acidi, alcalini o solventi come ad esempio l'acetone o il diluente	
Manutenzione e Pulizia Ruote	
Pulire le ruote con acqua tiepida e detersivo neutro	
Lubrificare periodicamente le ruote	
Conservazione e Smaltimento	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo MDR 2017/745	Classe I
Codice ISO DM332/99	12.06.09.009 / 12.06.09.106 / 12.06.09.109 / 12.06.09.130 / 12.06.09.133 / 12.06.09.142
Codice ISO DPCM 12.01.2017	12.06.09.003
Codice CND	Y120609
Codice RDM	2487867/R

Product Name  
Rollator Compacto.



Code	Colors
01.06.401A	Silver
01.06.401GO	Matte Gray
01.06.401B	Night blue

**Use**

The walkers and rollators are indicated in the walk rehabilitation for elderly or subjects with functional limitations of various kinds. They offer stability, ease of use and improve the sense of balance.  
This aid is designed for a single user, in indoor or outdoor environments according on models.

	(Seat Height Level 1) *	(Seat Height Level 2) *	(Seat Height Level 3) *	(Seat Height Level 4) *
Height of handle (cm)	77 to 90 (in 8 positions)	79 to 92 (in 8 positions)	81 to 94 (in 8 positions)	83.5 to 96.5 (in 8 positions)
Seat height (cm)	50	52	54	56,5
Seat dimensions (cm)	36 x 31	36 x 31	36 x 31	36 x 31
Length unfold (cm)	68,5	71	74	76,5
Width unfold (cm)	62	62	62	62
Distance between handrails (cm)	46	46	46	46
Wheel size (cm)	20 (8")	20 (8")	20 (8")	20 (8")
Shopping bag dimensions (cm)	33.5 x 21 x 15	33.5 x 21 x 15	33.5 x 21 x 15	33.5 x 21 x 15
Max. Load capacity of user (kg)	110 kg	110 kg	110 kg	110 kg
Net Weight (kg)	7,2	7,2	7,2	7,2
Gross Weight (kg)	7,8	7,8	7,8	7,8
Carton size	54.5 x 16 x 49.5	54.5 x 16 x 49.5	54.5 x 16 x 49.5	54.5 x 16 x 49.5

**Description**

- Foldable aluminum structure.
- Seat adjustable in 4 heights.
- Padded seat and back with lateral stick holder.
- Four wheels with a diameter of 20 cm (swiveling front wheels, two fixed rear wheels with brakes).
- Storage basket included.
- Ergonomic and height-adjustable knobs.
- Padded backrest adjustable in inclination.

Materials	
Frame	Aluminum
Knobs	Polypropylene
Wheels	PU
Seat	Padded PVC
Frame Maintenance and Cleaning	
Use a damp cloth with mild soap	
Do not use water and / or steam jet washing machines, abrasive objects and metal brushes	
Must not use acid, alkaline or solvent products such as acetone	
Wheels Maintenance and Cleaning	
Clean the wheels with warm water and neutral detergent	
Lubricate the wheels periodically	
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Class I
ISO Code DM332/99	12.06.09.009 / 12.06.09.106 / 12.06.09.109 / 12.06.09.130 / 12.06.09.133 / 12.06.09.142
ISO Code DPCM 12.01.2017	12.06.09.003
CND Code	Y120609

RDM Code	2487867/R
----------	-----------