

**Nome Prodotto**

Protesi Maxima Asimmetrica Piatta.



| Codice          | Lato     | Taglia | Taglia Francese | Taglia Europea |
|-----------------|----------|--------|-----------------|----------------|
| 06.01.6003PD.0  | Destro   | 0      | 65              | 50             |
| 06.01.6003PD.1  | Destro   | 1      | 70              | 55             |
| 06.01.6003PD.2  | Destro   | 2      | 75              | 60             |
| 06.01.6003PD.3  | Destro   | 3      | 80              | 65             |
| 06.01.6003PD.4  | Destro   | 4      | 85              | 70             |
| 06.01.6003PD.5  | Destro   | 5      | 90              | 75             |
| 06.01.6003PD.6  | Destro   | 6      | 95              | 80             |
| 06.01.6003PD.7  | Destro   | 7      | 100             | 85             |
| 06.01.6003PD.8  | Destro   | 8      | 105             | 90             |
| 06.01.6003PD.9  | Destro   | 9      | 110             | 95             |
| 06.01.6003PD.10 | Destro   | 10     | 115             | 100            |
| 06.01.6003PD.11 | Destro   | 11     | 120             | 105            |
| 06.01.6003PD.12 | Destro   | 12     | 125             | 110            |
| 06.01.6003PD.13 | Destro   | 13     | 130             | 115            |
| 06.01.6003PD.14 | Destro   | 14     | 135             | 120            |
| <hr/>           |          |        |                 |                |
| 06.01.6003PS.0  | Sinistro | 0      | 65              | 50             |
| 06.01.6003PS.1  | Sinistro | 1      | 70              | 55             |
| 06.01.6003PS.2  | Sinistro | 2      | 75              | 60             |
| 06.01.6003PS.3  | Sinistro | 3      | 80              | 65             |
| 06.01.6003PS.4  | Sinistro | 4      | 85              | 70             |
| 06.01.6003PS.5  | Sinistro | 5      | 90              | 75             |
| 06.01.6003PS.6  | Sinistro | 6      | 95              | 80             |
| 06.01.6003PS.7  | Sinistro | 7      | 100             | 85             |
| 06.01.6003PS.8  | Sinistro | 8      | 105             | 90             |
| 06.01.6003PS.9  | Sinistro | 9      | 110             | 95             |
| 06.01.6003PS.10 | Sinistro | 10     | 115             | 100            |
| 06.01.6003PS.11 | Sinistro | 11     | 120             | 105            |
| 06.01.6003PS.12 | Sinistro | 12     | 125             | 110            |
| 06.01.6003PS.13 | Sinistro | 13     | 130             | 115            |
| 06.01.6003PS.14 | Sinistro | 14     | 135             | 120            |

**Destinazione d'uso**

Le protesi mammarie esterne vengono utilizzate in soggetti che hanno subito una mastectomia. Grazie a questi dispositivi la morfologia del corpo resta invariata.

**Descrizione**

La PROTESI MAXIMA ASIMMETRICA PIATTA è composta da un gel di silicone avvolto da un sottile film di poliuretano; concepita per ripristinare una visione anatomico-estetica del corpo, specialmente per i pazienti che hanno subito un'asportazione del cavo ascellare. La protesi è anallergica, esteticamente naturale ed ergonomica grazie alla forma allungata. Sono disponibili sia la protesi destra che quella sinistra.

| Materiali  |  |
|--|--|
| Protesi  | Silicone   |
| Manutenzione e Pulizia   |  |
| La pulizia è possibile con un panno umido e un detergente neutro   |  |
| Si consiglia di utilizzare deodoranti, profumi o altri cosmetici prima di indossare le protesi                       |  |
| Conservazione e Smaltimento  |  |
| Si consiglia l'uso ad un solo utente   |  |
| Conservare al riparo da fiamme libere e fonti di calore  |  |
| Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata |  |
| Confezione e Nr. Articoli  |  |
| Scatola  | 1 Pezzo  |
| Garanzia   |  |
| Europea  | 24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili |
| Riferimenti del Prodotto   |  |
| Fabbricante  | Demarta Virginio S.r.l.  |
| Classe di Dispositivo MDR 2017/745   | Classe I   |
| Codice ISO   | 06.30.18.006   |
| Codice CND   | Y063018  |
| Codice RDM   | 1463825/R  |

**Product Name**

Asymmetrical Prosthesis Maxima with Flat Back.



| Code            | Side  | Size | French Size | European Size |
|-----------------|-------|------|-------------|---------------|
| 06.01.6003PD.0  | Right | 0    | 65          | 50            |
| 06.01.6003PD.1  | Right | 1    | 70          | 55            |
| 06.01.6003PD.2  | Right | 2    | 75          | 60            |
| 06.01.6003PD.3  | Right | 3    | 80          | 65            |
| 06.01.6003PD.4  | Right | 4    | 85          | 70            |
| 06.01.6003PD.5  | Right | 5    | 90          | 75            |
| 06.01.6003PD.6  | Right | 6    | 95          | 80            |
| 06.01.6003PD.7  | Right | 7    | 100         | 85            |
| 06.01.6003PD.8  | Right | 8    | 105         | 90            |
| 06.01.6003PD.9  | Right | 9    | 110         | 95            |
| 06.01.6003PD.10 | Right | 10   | 115         | 100           |
| 06.01.6003PD.11 | Right | 11   | 120         | 105           |
| 06.01.6003PD.12 | Right | 12   | 125         | 110           |
| 06.01.6003PD.13 | Right | 13   | 130         | 115           |
| 06.01.6003PD.14 | Right | 14   | 135         | 120           |
| 06.01.6003PS.0  | Left  | 0    | 65          | 50            |
| 06.01.6003PS.1  | Left  | 1    | 70          | 55            |
| 06.01.6003PS.2  | Left  | 2    | 75          | 60            |
| 06.01.6003PS.3  | Left  | 3    | 80          | 65            |
| 06.01.6003PS.4  | Left  | 4    | 85          | 70            |
| 06.01.6003PS.5  | Left  | 5    | 90          | 75            |
| 06.01.6003PS.6  | Left  | 6    | 95          | 80            |
| 06.01.6003PS.7  | Left  | 7    | 100         | 85            |
| 06.01.6003PS.8  | Left  | 8    | 105         | 90            |
| 06.01.6003PS.9  | Left  | 9    | 110         | 95            |
| 06.01.6003PS.10 | Left  | 10   | 115         | 100           |
| 06.01.6003PS.11 | Left  | 11   | 120         | 105           |
| 06.01.6003PS.12 | Left  | 12   | 125         | 110           |
| 06.01.6003PS.13 | Left  | 13   | 130         | 115           |
| 06.01.6003PS.14 | Left  | 14   | 135         | 120           |

**Use**

External breast prosthesis are used in subjects who have undergone mastectomy. Thanks to these devices the morphology of the body remains unchanged.

**Description**

The ASYMMETRICAL PROSTHESIS MAXIMA WITH FLAT BACK is composed of a silicone gel wrapped in a thin polyurethane film; designed to restore an anatomical-aesthetic vision of the body, particularly for those who went through a removal of the axillary cavity. The prosthesis is hypoallergenic, aesthetically natural and ergonomic thanks to the elongated shape. Both right and left prosthesis are available.

| Materials  |  |
|--|--|
| Prosthesis   | Silicone Gel   |
| Maintenance and Cleaning   |  |
| Cleaning is possible with a damp cloth and a mild soap   |  |
| Avoid direct contact with cosmetics (deodorant, perfumes...), please use them before putting the prosthesis on         |  |
| Storage and Disposal   |  |
| Single-user Product  |  |
| Keep away from heat, sparks, open flame or other sources of ignition   |  |
| Dispose of the device in compliance with current regulations on environmental protection and separate waste collection |  |
| Package and Nr. Items  |  |
| Box  | 1 Piece  |
| Warranty   |  |
| European   | 24 months for manufacturing defects with the exception of consumable parts |
| Product References   |  |
| Manufacturer   | Demarta Virginio S.r.l.  |
| Directive MDR 2017/745   | Class I  |
| ISO Code   | 06.30.18.006   |
| CND Code   | Y063018  |
| RDM Code   | 1463825/R  |